

# SKUP evaluations CRP, PT-INR, HBA1c and Haemoglobin

(SKandinavisk Uprøvning af laboratorieudstyr til Primærsektor)

Esther A Jensen, SKUP, Hillerød Hospital. Per Grinsted, DAK-E [www.skup.dk](http://www.skup.dk) [www.skup.nu](http://www.skup.nu)

Komponenter afprøvet:

Allergentest, CRP, Glucose, HBA1c, Hæmoglobin, Hæmatologi, Leukocytter + 3-parts diff., Mononucleose, PT-INR, SR, Strep A, U-albumin, U-hCG

Kvalitetsmål



Total Error

Bias

Imprecision

Glucose

20 %

< 3 %

< 4 %

Hb

5 %

< 2 %

< 3 %

HBA1C

10 %

< 4 %

< 4 %

PT-INR

20 %

< 6 %

< 5 %

CRP

23 %

<10 %

<10 %

Leukocyt

16 %

-

-

U-HCG

< 5 eller > 40 IU/L

< 5 eller > 40 IU/L



HBA1c

Hospital

Lægehuse



HB

Hospital

Lægehuse



Afinion

SKUP/2008/65





Biotest

SKUP/2001/17




DCA 2000

SKUP/1999/4




ABX Micros

SKUP/2002/23\*



Hemo\_Control

SKUP/2004/29



Chempaq

SKUP/2006/47



PT-INR

Hospital

Lægehuse



CRP

Hospital

Lægehuse



CoaguCkek S

SKUP/2000/7




QuickRead

SKUP/2001/12




Thrombotrack

SKUP/2000/8




ABX Micros

SKUP/2002/23\*





Pro Time

SKUP/2000/11





HemoChron

SKUP/2004/33




i-CHROMA

SKUP/2007/61



CoaguCkek SX

SKUP/2007/55



Simple Simon PT

SKUP/2007/57\*

